BIOPRESS™ CULTURE PLATES

Flexible bottomed culture plate for use with the Flexcell Compression system.

- Each well contains a stationary platen and piston between which 3-D tissue samples or cell-seeded constructs are deformed using the FX-5000 Compression system.
- Ring foam holders keep samples confined to central region of each well.
- All plates come pre-sterilized.
- Available individually or in cases of 40 plates.



	ORDERING INFORMATION (please call your discipation for prices)	
Catalog Number	Product/Item	
BF-3000C	BioPress Compression Culture Plate	

世联博研北京科技公司-美国FKANI在中国大陆独家授权代理商