


**Disposable Rodent Caging System from Innovive** is the top-of-the-line research-ready rodent housing solution.

All caging components for Innovive Disposable Rodent Caging are manufactured from post-consumer PETE plastic using solar energy. Cages arrive at your facility double bagged and irradiated, eliminating the need for any washing or sterilizing. They are packaged as 4 bags of 25 cages for a total of 100 per box.

**Good for the environment. Good for your science.**  
**BPA-FREE**

 Made in USA

Innovive Disposable Rodent Caging utilizes a **universal cage bottom** for its line of mouse caging solutions. The cage bottom offers 81 in<sup>2</sup> of floor space and a minimum 5" of headroom. All cage lids form a secure air-tight seal when snapped onto the cage bottom.

Designed for compact storage, Innovive caging offers dense stacking: a stack of 25 cage bottoms is only 13" tall; 25 changeouts fit easily into a standard changing hood.

Cage bottoms and lids are sold as one unit. Cage bottoms can also be purchased separately.

 green product



### Innocage® Static Mouse

14.7"D x 9.2" W x 5.5" H | 8.8 oz cage | 1.9 oz food tray

**Watering:** Aquavive® Mouse bottle

The Innocage® Static Mouse consists of a universal mouse cage bottom and a static cage lid. Add an Innovive food tray, water bottle and cage card holder to build a complete cage set-up.

No IVC rack needed. Ideal for short-term studies or emergency preparedness.

- Large filter areas create a standard micro-isolator environment.
- External access to water bottle.
- 300 ml water bottle (Aquavive®) allows for 2 week- changes
- High capacity, extended - use feeder (sold separately)

 green product



### Innocage® Transport Mouse

14.7"D x 9.2" W x 5.5" H | 8.8 oz cage

**Watering:** gels

The Innocage® Transport is a great solution for moving animals a long or short distance: it keeps mice comfortable, secure and well ventilated throughout their journey. The cage consists of a universal cage bottom and a lid designed specifically for rodent transportation.

- Low-cost transportation solution.
- Secure snap-top lid.
- Disposable cages eliminate tracking inventory of expensive washable cages from your operations.
- Optional feeder and card holder (sold separately)

**NOTE:** All Innovive cages packaged 100 units per box (4 bags of 25 lids). Cage purchase includes lids and cage bottoms (bottle, feeder and card holders sold separately). Min. cage order is 100 units.

