


HT BIOFLEX™ 细胞培养板

HT BIOFLEX™ CULTURE PLATES

24 孔高质量细胞培养板

- 24孔圆形,总生长面积37.47 cm² (1.56 cm²/孔).孔板底部是由弹性硅胶膜制作而成,柔韧性好,透明度高.
- 孔板底部的硅胶膜厚度仅 0.254 mm,可以使用标准正立式显微镜观察细胞.
- 可以施加最高15% 的等量双轴向应力.
- 五种不同包被的培养表面: Amino, Collagen (Type I or IV), Elastin, ProNectin (RGD), Laminin (YIGSR).
- 两种颜色可供选择:白色和黑色.



定货信息(请与销售商或公司联系询问价格)

产品编号	英文名称
细胞培养板颜色: 白色	
HTPW-3001U	White HT BioFlex Culture Plate — Untreated
HTPW-3001A	White HT BioFlex Culture Plate — Amino
HTPW-3001C	White HT BioFlex Culture Plate — Collagen Type I
HTPW-3001C(IV)	White HT BioFlex Culture Plate — Collagen Type IV
HTPW-3001E	White HT BioFlex Culture Plate — Elastin
HTPW-3001P	White HT BioFlex Culture Plate — ProNectin
HTPW-3001L	White HT BioFlex Culture Plate — Laminin
细胞培养板颜色: 黑色	
HTPB-3001U	Black HT BioFlex Culture Plate — Untreated
HTPB-3001A	Black HT BioFlex Culture Plate — Amino
HTPB-3001C	Black HT BioFlex Culture Plate — Collagen Type I
HTPB-3001C(IV)	Black HT BioFlex Culture Plate — Collagen Type IV
HTPB-3001E	Black HT BioFlex Culture Plate — Elastin
HTPB-3001P	Black HT BioFlex Culture Plate — ProNectin
HTPB-3001L	Black HT BioFlex Culture Plate — Laminin