

BIOPRESS™ CULTURE PLATES

Flexible bottomed culture plate for use with the Flexcell Compression system.

- Each well contains a stationary platen and piston between which 3-D tissue samples or cell-seeded constructs are deformed using the FX-5000 Compression system .
- Ring foam holders keep samples confined to central region of each well.
- All plates come pre-sterilized.
- Available individually or in cases of 40 plates.



ORDERING INFORMATION (please call your distributor for prices)

Catalog Number	Product/Item
BF-3000C	BioPress Compression Culture Plate